

APPROVED	O.G. FIG.	
BY	CLASS	SUBCLASS
CRAFTSMAN		

Instructions

- 1) Do not eat or drink 1 hour prior to collection of the DNA sample.
- 2) Wash and dry hands thoroughly.
- 3) Remove the contents of the collection kit.
- 4) Remove the storage card and desiccant bag from the resealable pouch.
 CAUTION! Do not touch the specimen area (circles) of the storage card.
 Do not touch the collection tip of the swab.
- 5) Complete the information section on the storage card.
- 6) Open the storage card to expose the specimen circles.
- 7) Rub the inside of the cheek and gums briskly with the swab.
- 8) Blot a specimen circle of the storage card.
- 9) Repeat steps 7 and 8 for each specimen circle on the storage card.
- 10) Allow the storage card to air dry for a minimum of 3 hours.
- 11) Return the storage card and desiccant bag to the resealable pouch and
- 12) Place the tamper-resistant seal over the resealable portion of the storage
- 13) Return storage pouch to the permanent storage record folder.
- 14) Store the folder at room temperature in a safe, dry location such as a drawer or safety deposit box.

Fig. 1

APPROVED	O.G. FIG.	
BY	CLASS	SUBCLASS
DRAFTSMAN		

DNA™ COLLECTION AND STORAGE KIT

Identification Information

Name: _____

Date of Birth: _____ - _____ - _____

Social Security Number: _____ - _____ - _____

Date of Collection: _____ - _____ - _____

Signature: _____

UPON MY DEMISE I BEQUEATH THIS SAMPLE TO:

Name: _____

Relationship: _____

Witness: _____

Witness: _____

Optional

<i>SUBSCRIBED AND SWORN TO BEFORE ME</i> this	
_____ day of _____,	_____
<i>Notarial Seal</i>	_____
	Notary Public in and for
	_____ County, _____
	My Commission Expires:

Fig. 2